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HYALONE®

Hyaluronic acid sodium salt solution for intra-articular injection

Pre-filled syringe 60 mg/4 ml

EN | Product description

HYALONE® is a sterile, viscoelastic solution manufactured with hyaluronic acid sodium salt, obtained by bacterial fermentation from a fraction of high molecular weight. Hyaluronic acid, a polysaccharide of the glycosaminoglycan family, is naturally present in many human tissues such as cartilage and synovial fluid; it is continuously secreted into the joint space and represents a major component of the synovial fluid, to which it provides its characteristic viscosity and elasticity. Such properties are fundamental for the lubricating and shock absorbing functions exerted by the fluid in normal joints to protect cartilage and soft tissues against mechanical injuries.

In traumatic and degenerative joint disorders, an insufficient amount of hyaluronic acid and a loss of viscosity occur in synovial fluid, resulting in an impairment of joint function and in a painful symptomatology. Extensive data in the literature indicate that intra-articular administration of hyaluronic acid is capable to restore the visco-elastic properties of the synovial fluid, with alleviation of pain and improvement of joint mobility. HYALONE® is a treatment for temporary synovial fluid replacement in patients affected by degenerative or mechanical arthropathy of the knee and hip, which causes an alteration of the functional performances of the synovial liquid. Intra-articular injection of HYALONE® reduces pain symptoms and improves knee and hip joint functionality, up to 6 months.

Composition

Principal component: Hyaluronic acid sodium salt 1.5%
Other components: Sodium chloride, Disodium hydrogen phosphate dodecahydrate, Sodium dihydrogen phosphate dihydrate, Water for injections.

Indications

HYALONE® is indicated for the treatment of pain and improvement of joint functionality in patients affected by degenerative or mechanical arthropathy of the hip and knee.

The product is indicated in the treatment of persistent pain after initial failure of analgesics or in case of failure or intolerance to non-steroidal anti-inflammatory drugs.

Intended use

HYALONE® is a temporary synovial fluid replacement to be administered by intra-articular injection.

Dosage and Administration

HYALONE® is intended for intra-articular injection only. Product administration should be performed exclusively by licensed physicians (e.g. orthopaedic surgeon, rheumatologist, physiatrist radiologist, sports doctors, etc.). Injection must be strictly intra-articular. Intra-articular injection has to be performed according to the usual standard technique, using precise, anatomical localization. In the hip, ultrasound guidance or radioscopy guidance is recommended.

Remove any joint effusion, if present, before the administration. All the rules regarding the aseptics and the injection technique shall be followed.

Connect a suitable sterile needle (for example 18 or 20 G) to the syringe and inject HYALONE® in the affected joint.

Further treatments after the first application may be needed to maintain the benefit of the treatment over time, depending from the individual patient needs.

Clinical data show that the effectiveness of HYALONE® persists up to 6 months. A further injection, if clinically requested, should be administered after this period of time.

Contraindications

Do not administer to patients with ascertained individual hypersensitivity to the product components.

Do not administer in cases of infections or skin diseases in the area of the injection site.

Do not administer to patients with active synovitis.

Use in specific populations

The safety and efficacy of HYALONE® in pregnant women, lactating women who are breast feeding, or in subjects less than 18 years of age have not been established and therefore its use is contraindicated in these patient populations.

Warnings and Precautions

Do not use in case of package damage.

Do not use the product after the expiry date reported on the package.

The expiry date refers to the product kept in its original package at a temperature not exceeding 25° C.

The syringe is for single use, that means it should be used once only for a single patient. Inject the contents in one joint only.

For the first 24 hours after injection, the patient is allowed to continue all routine activities of daily living but is recommended not to overuse the treated joint.

The assembled syringe must be discarded immediately after use, regardless of whether or not the solution has been completely

administered.

If this product is reprocessed and/or reused, Fidia Farmaceutici cannot guarantee performance, functionality, material structure, or cleanliness or sterility of the product. Reuse could lead to illness, infection and/or serious injury to the patient or user.

After use, dispose according to applicable national practice.

Keep out of reach of children.

Undesirable effects

Local pain, swelling, heat and redness may occur sporadically at the injection site. Such symptoms are generally mild and transient. Following intra-articular injection, application of an ice pack onto the treated joint for five to ten minutes will reduce the incidence of these events. Local or systemic allergic reactions may occur in subjects with hypersensitivity to the product components.

More marked inflammatory reactions, sometimes with calcium pyrophosphate crystals, have been occasionally reported in association with intra-articular injections of hyaluronate.

As for any intra-articular treatment, septic arthritis may rarely occur when general precautions for injections are not observed or the site of injection is not aseptic.

Interactions

Do not use concomitantly with disinfectants containing quaternary ammonium salts, because hyaluronic acid can precipitate in their presence.

In order to prevent any possible interactions, avoid the contemporary administration of HYALONE® with other intra-articular products.

Storage

Store at a temperature not exceeding 25° C.

How supplied

Box containing one pre-filled syringe with 60 mg/4 ml solution of hyaluronic acid sodium salt. The solution is clear and transparent. The eventual presence of a bubble does not compromise the quality of the product.

The syringe, sterilised using steam, is pre-assembled with a backstop and is sealed in a blister.

Available packaging:

- Box containing 1 pre-filled syringe
 - Box containing 50 pre-filled syringes (for healthcare facilities only)
 - Box containing 100 pre-filled syringes (for healthcare facilities only)
- Needles are not included in the package.

Manufacturer FIDIA Farmaceutici S.p.A.

Via Ponte della Fabbrica 3/A
35031 Abano Terme (PD) Italy

Date of the latest revision of the instructions for use: January 2025

Consult instructions for use

Use by date

Batch code

Do not reuse

Sterilised using steam

Single sterile barrier system

Single sterile barrier system with protective packaging outside

Do not use if package is damaged

Temperature upper limit

Manufacturer

CE Mark and Notified Body number