#### STRUMENTARIO CHIRURGICO

#### ITALIANO

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#### SURGERY TOOLS

#### ENGLISH

The product can be used by specialized medical staff only who have been taught how to handle the tools. Improper use or little or inappropriate maintenance can rapidly damage the tools. The first and anytime after use, the product must be thoroughly cleaned with warm water, disinfected, rinsed (if possible with distilled water) and carefully dried. It must be sterilized before the next use. The product can be cleaned in special ultrasonic tanks and disinfected in tool washers with a temperature of maximum 95°C (203°F). After cleaning, and before sterilizing, it is advisable to handle the tools with physiologically safe oil, above all as regards the tips, the connectors, the terminals and all movable parts. Moreover, pay attention that the product does not come into contact with acids or other aggressive disinfectants which might corrode it.

The advised sterilization method is through autoclave with steam, with the temperature of the cycle not exceeding 135°C (275°F) for maximum 15 minutes, in order to avoid damaging the product. For warm air sterilization, it is advisable to use a temperature between 180° and 200°C. The products can be sterilized in tools cleaning machines, too, up to a temperature of 121°C. The tools showing corrosion signs must be immediately thrown away.

In guarantee it will be possible to freely repair and/or replace all those parts presenting well proved production defects, except for the expenses for labour, transfer, transport, packaging etc. Damages due to improper use of the product cannot be covered by the guarantee.

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